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- Castell, W.P. et al. Cholesterol and Other Lipids in Coronary Heart Disease Circulation 55:757
- Stein, E.A. Myers, G.L. Lipids, Lipoproteins and Apolipoproteins in Tietz Textbook of Clinical Chemistry, Second Edition, Burtis, C. Ashwood, A. (Eds.) CV Mosby Co. 23:1002-83 (1994).
- Warmick, R. Nguyen, T. Albens, A.A. Comparison of Improved Precipitation Methods for Quantification of High Density Lipoprotein Cholesterol - Clin. Chem. 1985: 31:217-22.
- Warmick, R. Wood, P. National Cholesterol Education Program Recommendations for Measurement of High Density Lipoprotein Cholesterol, Executive Summary, Clinical Chemistry, Vol. 41, N° 10, 1996.
- National Committee for Clinical Laboratory Standards. National Evaluation Protocols for Interference Testing, Evaluation Protocol Number 7, Vol. 4 N° 8, June 1984
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, Third Edition, AACC Press (1990)
- National Institute of Health Publication N° 01-3670, May 2001.
- Arquivos da EBRAM.

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GENERAL	
Measurement Mode:	ABSORB
Reaction Mode:	R-S-SR1
Calibration Mode:	Calibration (Single AVG)
Reagent Blank:	Reagent/Dil
Cleaner:	Before
Wavelength:	600 nm
Decimal Position:	0
Unit:	mg/dL
ANALYSIS	
Dilution Name:	STD
Factor:	Main STD
Time:	Sample Dil Name: H2O
Post Dil Factor: 2.0	Conc. Factor: No
Sample:	Circle: 1
Volume:	Dil: 10 ul
Reagent:	Circle: 1
Volume:	200 ul
Start Reagent 1:	Circle: 12
Volume:	Dil: 5.0 ul
Start Reagent 2:	Circle:
Volume:	Dil:
CALCULATION	
Sample Limit:	NO
Point:	
Reaction Direction:	Increase
Check:	On
Conversion Factor:	1.0000
Offset:	0.0000
Test Range Low:	1.06
High:	100
Norm Range Low:	30
High:	85
Number of Slats:	1
Calculation Step A:	Endpoint
Readings First:	11
Last:	24
Reaction Limit:	
Point:	
Calib. Interval:	On Request
CALIBRATION	
Reagent Blank:	
Reagent Range:	
Reagent Range Low: NO	Blank Range Low: NO
High: NO	High: NO
Factor:	
Calibrator Pos:	*
STD1:	**
STD2:	
STD3:	
STD4:	
STD5:	
STD6:	
STD7:	
STD8:	
Calc. Model:	
Correction STD:	
Replicate:	Duplicate
Deviation:	5%
CONTROL	
CS1 - Pos (?) Assay (***) Low (***) High (***)	
CS2 - Pos (?) Assay (***) Low (***) High (***)	
CS3 - Pos (?) Assay (***) Low (***) High (***)	
(*) colocar a posição correspondente da rack CAUCS	
(**) colocar o valor correspondente do calibrador	
(***) colocar o valor correspondente do soro controle	

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Test Name: HDL Colestamol	Test: HDL		
Test Bar Code:			
Test Type: endpoint	Curve Type: Blank Linear		
Units: mg/dl	Nº of Decimal Places: 0		
Primary Wavelength: 600	Secondary Wavelength: none		
Read Time Interval: 20	Sample Blank: Yes		
Factor:			
Calibration Interval: (*)			
Normalization Interval:			
Nº of Calibrators: 2	Nº of Replicates: 2		
Low Blank A Limit: -0.010	High Blank A Limit: 0.350		
Low A Limit: -0.050	High A Limit: 2.000		
Low Normal: 30	High Normal: 85		
Linearity Limit: 100	Curve S.D. Limit: 5.000		
Test Name: HDL Colestamol	Test: HDL		
Test Bar Code:			
Sample Volume: 3 ul	Sample Diluent:		
Reagent Dilution Ratio: 1	Predilution Ratio: 1		
Reagent Dilution:			
Reagent Volume	Bar Code	Diluent Volume	Lag Time
Reagent 1	250	YH1A	290 sec
Reagent 2	60	YH2A	290 sec
Reagent 3			
Reagent 4			
Controls:			
(*) Introduzido pelo operador (1 - 999)			

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Name:	HDL Colestamol
Abbr. Name:	HDL COL
Mode:	Endpoint
Wavelength:	546 nm
Units:	mg/dL
Decimals:	0
Low Conc:	1.06
High Conc:	100
Calibrator Name:	*
Repeat:	2
Number:	1
Concentration:	**
Interval:	*
Cut off:	NO
Prozone Check:	NO
Ref Male Low:	30
Ref Male High:	70
Ref Female Low:	30
Ref Female High:	85
Ref Pad Low:	*
Ref Pad High:	*
Control 1:	*
Control 2:	*
Control 3:	*
Correlat Factor:	1.000
Correlat Offset:	0.000
DUAL MODE	
Name:	HDL Colestamol
Sample Blank:	NO
R1 Bottle:	
Normal Volume:	300 ul
Reagent Volume:	301 ul
SAMPLE	
Normal Volume:	3 ul
Reagent Volume:	2 ul
R2 Bottle:	
Normal Volume:	100 ul
Reagent Volume:	100 ul
Predilution:	NO
Slope Blank:	
Delay Min Time:	
Linearity Limit:	
Point One, Two:	
Incubation Time:	6.5 min
Low Absorbance:	-0.100
High Absorbance:	3.000
R ABS L Limit:	-0.100
R ABS H Limit:	3.000
Substr. Depletion:	
Reagent Blank:	
R ABS Deviation:	YES #
Cal Low Limit:	*
Cal High Limit:	*
Factor:	
MONO MODE	
Name:	
Sample Blank:	
R1 Bottle:	
Normal Volume:	
Reagent Volume:	
SAMPLE	
Normal Volume:	
Reagent Volume:	
Delay Min Time:	
Linearity Limit:	
Predilution:	
Incubation Time:	
Point One, Two:	
Low Absorbance:	
High Absorbance:	
R ABS L Limit:	
R ABS H Limit:	
Substr. Depletion:	
Reagent Blank:	
R ABS Deviation:	
Reagent Blank:	
Cal Low Limit:	
Cal High Limit:	
Factor:	#
(*) Dados colocados pelo usuário	
(**) Colocar o valor correspondente do calibrador	
(***)	
(#) Dados calculados pelo analisador	

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Name:	HDL - Chol
Immunoassay:	NO
Chemistry Type:	Endpoint
BioType:	
Innova Chemistry:	NO
% Sample Volume:	6 (3.0 ul)
Wavelength:	560
Biochromatic Chemistry:	NO
Biochromatic Factor:	
K1:	
K2:	
Bio Limit 1:	
Bio Limit 2:	
Resolution Limit:	
Delay Time:	5.00
Blank Type:	NO Blank
Incubation:	
% Reagent Volume:	50 (250 ul)
2nd Reagent:	YES
2nd Reagent Volume:	18 (90 ul)
A2 Delay:	3.00
Units:	mg/dL
Unit Factor:	1.000
Decimal Point:	0
RBL Low:	0
RBL High:	0.50
Range Low:	1.06
Range High:	100
Calibration Factor:	(*)
Reagent Rate:	
Standard Value:	(**)
Normal Low:	30
Normal High:	85
Slope:	1.0000
Intercept:	0.0000
C1*10E-6:	
C2*10E-6:	
D1*10E-6:	
Delta R:	
Linear Factor:	
Final Limit:	
Endpoint Limit:	0.5000
(*) Colocar 1 como valor inicial. O valor do CAL FACT é determinado pelo ensaio de calibração	
(**) Colocar o valor correspondente do calibrador	

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Measure Unit:	mg/dL	Decimals:	0	Reaction Type:	Endpoint
Volume (ul):	Serum: 3	Plasma:		Units:	
Abs Range (m Abs):	Min: 0	Max: 2.000		Reagent 2: 100	
Linearity Limit:	100	Reagent Blank:	Yes		
Contaminating:	Yes	Differential:	No		
Filter 1:	546	Filter 2:	none		
Time (sec):	Mix 1: 2.00	Incubation 1: 300	Lag Phase: 2		
	Mix 2: 2.00	Incubation 2: 180	Measure: 1		
Measurement Type:	Standard	Factor:	0.00		
Normal Range					
Age (yrs)		Homens		Mulheres	
Below 10:	Min.	Max.	Min.	Max.	
From 10 to 80:	30	70	30	85	
Over:					

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