QUIMIPROT – Proteína Total

Biureto

ANVISA: 10159820240

Finalidade

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ANVISA: 10159820240

Purpose

Colorimetric reaction for the quantitative determination of total protein in human serum, ascitic, pleural and synovial fluid samples. For "in vitro" diagnostic use only.

Principle

The protein present in the sample reacts with copper (II) ions in an alkaline medium, creating a colored complex that can be quantified by spectrophotometry (550 nm). The intensity of the violet color is proportional to the total protein present, when compared to a solution with a known protein concentration.

Methodology

Biuret

Clinical significance

The level of total protein in serum is basically a reflection of hepatic synthesis or protein loss due to kidney disease, severe malnutrition, severe burns and hemodilution.

Increased levels are found in dehydration, multiple myeloma, chronic liver disease, macroglobulinemia, cryoglobulinemia, lupus erythematosus, rheumatoid arthritis, sarcoidosis, chronic infections, lymphogranuloma and sub-acute bacterial endocarditis.

Decreased levels are found in kidney disease, nephrosis, hyperhydration, severe malnutrition, severe burns, malabsorption syndrome, calcium and vitamin D deficiency, and terminal liver disease.

Ascitic fluid: values below 2.5 g/dL are found in cirrhosis and heart failure, above 3.0 g/dL are found in carcinomatosis, chylous ascites, pancreatitis.

Pleural fluid: values below 2.5 g/dL are found in cirrhosis, heart failure and nephrotic syndrome and above 3.0 g/dL are found in neoplasms, infections, pancreatitis, collagenosis, embolism and chylothorax.

Synovial fluid: high values can occur in inflammatory joint processes.

Reagents

Single reagent ready for use. Store at 2 – 25°C and keep away from light. Contains: 600 mM sodium hydroxide, 12 mM cupric sulphate, 32 mM sodium potassium tartrate, 30 mM potassium iodide.

Standard (code: 3006): Store at 2 – 8°C. Aqueous solution containing total protein standard concentration traceable to Standard Reference Material ERM-DA470k/IFCC. Check the concentration of the standard on the vial label.

Unopened reagents are stable until the expiration date printed on the product label and on board (in a refrigerated compartment of the analyzer) have a stability of approximately 20 days. During handling, the reagents are subject to chemical and microbial contamination which can lead to a reduction in stability.

Precautions and precautions required

This reagent should only be used for in vitro diagnostics. The reagent contains sodium hydroxide which is corrosive and can cause burns. Do not pipette by mouth. Avoid contact with skin and clothing. In case of contact with eyes, rinse with plenty of water and seek medical attention.

The standard contains sodium azide as a preservative (0.05%). This component can react with copper and lead and can become an explosive metal. When disposing of it, add a large quantity of water.

The temperature of the working environment and the reaction time must be monitored in order to obtain correct results. Do not use if the absorbance of the blank exceeds 0.150 (converted to 1.0 cm of optical space) when measured at 550 nm and if it is difficult to achieve the values established for the fresh control serum, which indicates contamination of the reagent. The presence of a precipitate indicates deterioration of the reagent and should not be used.

Material required not supplied

1. Water bath or analyzer capable of maintaining a temperature of 37°C and capable of measuring absorbance at 550 nm (540 – 560nm).
2. Pipettes for measuring samples and reagent.
3. Distilled/deionized water.
4. Analyzer consumables when used.
5. Control sera and calibrators.
6. Stopwatch.

Sample

Serum – the protein in the sample is stable for 3 days if refrigerated at 2 – 8°C or 7 days at -20°C. Ascitic, pleural and synovial fluids – centrifuge for 10 minutes at 3000 rpm and use the supernatant to proceed with the test.

Do not use urine or cerebrospinal fluid samples, due to the methodology used.

The use of a tourniquet for more than 3 minutes causes an increase in the albumin value.

All samples and controls are considered potentially infectious, so we suggest handling them following the established Biosafety standards

Patient preparation

A 6-hour fast and collection in the morning is recommended. However, this may be modified according to medical advice.

Interferences

- Bilirubin ≥ 20 mg/dL, Hemoglobin ≥ 100 mg/dL, Triglycerides ≥ 196 mg/dL may interfere with the result.
- Lipemia, gross hemolysis or plasma expanders (Dextran, PVP and Hemacel) can cause falsely high results. All anticoagulants interfere with the results.
- Some drugs and substances affect the protein concentration, we suggest consulting Young et al.

System parameters:

Temperature	37°C
Wavelength	550 nm (540 – 560 nm).
Reaction Type	End-point Colorimetric
Direction	Increasing
Sample x Reagent Ratio	1:100
Sample Vol.	10µL
Reagent Vol.	1.0 mL
Incubation Time	10 minutes

Calibration

Use Quimicalib Ebram Code 7023/12023 or the standard that comes with the kit, for kits with a standard. The concentration of total protein in Quimicalib is traceable to the NIST reference material SRM927d and in the standard to ERM-DA470k/IFCC.

Automated procedure

Application in the automated system: see manual for use of the equipment and instructions for use of the reagent.

Application in the semi-automatic system: proceed as shown below in the manual procedure up to item 2 (incubation), then use the equipment for reading, following the specific analytical protocol based on the System Parameters item.

Manual procedure

1. Separate 3 test tubes and carry out the procedures as follows:

	White	Standard	Sample/S.C.
Distilled water	10µL	-	-
Calibrator		10µL	
Sample	-	-	10µL
Reagent	1.0 mL	1.0 mL	1.0 mL

2. Homogenize the tubes and leave in a water bath (WB) at 37°C for 10 minutes. The water level in the water bath must be higher than the reagent level in the test tubes.
3. Zero the device with the reagent blank at 550 nm (540 – 560nm), take the readings and record the absorbance of the standard, sample and control serum (SC).

Note: the final color is stable for 60 minutes at room temperature.

* Strongly lipemic sera require a sample blank. Add 0.010 mL (10 µL) of sample to 1.0 mL of saline solution, homogenize and read the absorbance against water. Subtract this value from the patient's absorbance to obtain the corrected reading. Note: This procedure is suggested for spectrophotometers which require a minimum volume of 1.0 mL and can be adjusted proportionally without affecting the performance of the test. Please note that sample volumes smaller than 10 µL increase measurement inaccuracy in manual applications.

Calculations

(Abs.=Absorbance)
(Conc. = Concentration)

$$\frac{\text{Total Protein Sample (g/dL)}}{\text{Sample (g/dL)}} = \frac{\text{Abs. Sample}}{\text{Abs. Standard}} \times \text{Conc. of standard (g/dL)}$$

Formula for calculating Globulin:

Globulin (g/dL) = Total Protein (g/dL) - Albumin (g/dL)

Example:

Absorbance of sample = 0.350
Absorbance of standard = 0.400
Concentration of standard = 8.0 g/dL

$$\text{Total protein} = \frac{0.350}{0.400} \times 8.0$$

Total protein of sample = 7.0 g/dL

Linearity

When performed as recommended, the test is linear up to 15 g/dL. Samples above 15 g/dL should be diluted with saline solution to between 0.1 and 15 g/dL and the re-sults multiplied by the dilution factor.

Quality control

Each laboratory should maintain an internal quality program that defines objectives, procedures, standards, tolerance limits and corrective actions. A defined system for monitoring the analytical variation of the measurement system should also be maintained. We recommend using the Quimicon-trol Normal and Quimicontrol Abnormal Ebram control serums Codes 7024/12024 and 7031/12031.

Expected values

Serum= 6.5 - 8.0 g/dL
Ascitic fluid - transudate = < 2.5 g/dL
Ascitic fluid - Exudate= > 3.0 g/dL
Pleural fluid - Transudate = < 2.5 g/dL
Pleural fluid - Exudate = > 3.0 g/dL
Synovial fluid = 2.5 - 3.0 g/dL
These values are given solely as a guide. It is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI): g/L Total Protein (g/dL) x 10 = Total Protein (g/L)

Comparative studies

Studies carried out between this procedure and a similar methodology produced the following sta-tistical results:

Number of samples	40
Range of results	6.24 - 8.11
Correlation Coefficient	0.990
Slope	0.999
Intercept	0.19 (g/dL)

Precision

Precision studies were carried out with two levels (normal and pathological) with each sample proces-sed 10 times and the following statistical data were found:

N=10	Level 1	Level 2
Mean (g/dL)	6.24	8.11
D.P. (g/dL)	0.06	0.11
C.V. (%)	0.9	1.4

Accuracy

Samples were processed for 10 consecutive days, once a day and in duplicate. The following statistics were found:

N=10	Level 1	Level 2
Mean (g/dL)	6.24	8.11
D.P. (g/dL)	0.10	0.13
C.V. (%)	1.6	1.6

Methodological Sensitivity

0,1 g/dL

Specificity

As the samples were randomly selected from outpatients and hospitalized patients, it can be inferred that the method has adequate methodological specificity.

Observations

- Proper cleaning and drying of the material used are fundamental factors for the stability of the reagents and obtaining correct results.
- The water used in the laboratory must be of the appropriate quality for each application. Thus, to prepare reagents and use in measurements, it must have resistivity ≥1 mega ohm or conductivity ≤ 1 microsiemens and silicate concentration < 0.1 mg/L (type II water). For rinsing the glassware, the water can be type III, with resistivity ≥ 0.1 megaohms or conductivity ≤ 10 micro-siemens. Use type II water for the final rinse.

Presentation

General Biochemistry Line: 1 x 200mL + 1 x 10mL
Bulk Line: 1 x 500mL

Bibliographic reference

- Tietz, Textbook of Clinical Chemistry, WB Saunders, Philadelphia (1986)
- Riegler, E. Anal. Chem. 53:242 (1914)
- Kingsley, G.R. J. Biol. Chem. 131:197 (1939)
- Kingsley, G.R. J. lab Clin. Med. 27:640 (1942)
- Weichselbaum, T. Amer. J. Clin. Path. 16:40 (1948)
- Gomall, A. et al. J. Biol. Chem. 177:752 (1949)
- Henry, R.J. et. Al. Clinical Chemistry: Principals and Technics, Harper & Row, New York, p. 415 (1974)
- Young, D.S. et al, Clin. Chem. 21:1D (1975)
- Tietz, N.W. Fundamentals of Clinical Chemistry. W.B. Saunders, Philadelphia p. 299 (1976).
- Miller,O.,Gonçalves,R.R.,Laboratório para o Clínico, 8 ed.,Atheneu,(1998).
- EBRAM Archives.

Revision: March 2026.

Universal symbols used on packaging in vitro diagnostics			
	See instructions for use	 Reagent	 Manufactured by
	Contents are sufficient for <n> tests	 Expiry date (last day of the month)	 Batch number
	Temperature limit (keep at)	 Product for diagnostics <i>in vitro</i>	 Catalogue number

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